

CHAMBER		
PENNING GAUGES		
MAIN CONSOLE		
GUN MOVEMENT		
KEY BOARD		
COLUMN SHELL		
SCINTILLATOR CAP		

NO. 1 OF 2 PAGES

PAGE TWO OF ANASPEC'S QEMSEM ISI S.E.M. SERVICE SCHEDULE

VACUUM VALVE		
ROTARY PUMP		
ANTI- SUCK BACK		
GUN VALVE		
FARADAY CUP		
NOISY FANS		
BACKLASH (ABT55)		
BLOWN BULBS		
CRT. BRIGHTNESS		
STAGE LUBE		
COLUMN VALVES		
CLEAN AIR LOCK		

OPERATIONAL CHECKS

Redo all pre-service checks except step 12.

	GOOD	COMMENTS
THERMAL SWITCH		
VACUUM		
FILAMENT STABILITY		
GUN ALLIGNMENT		
STAGE RESOLUTION		
BACKSCATTER DETC.		
RESOLUTION		
LN2 LEVEL SENSE		
E.D.S. RESOLUTION		

BSD Au / FARADAY CUP		232 and 5
----------------------	--	-----------

GENERAL COMMENTS

SERVICE TECH

CUSTOMER SIGNATURE

.....

.....